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## Complete denture prosthetics with addition curing silicones

**General information:**

For each field of dental impression, Bisico offers a varied range of materials specifically adapted to your needs. Also in complete denture prosthetics, addition-curing silicones prove to be more efficient than other systems. Generally, they are kneaded at a 1:1 ratio by hand or are mixed with a spatula on the pad. Avoid contaminations between the components by mixing up trays and lids. Recommendations for functional impressions with Bisico.

**Function:**

The essential element of a functional impression is a precise functional tray adjusted to the specific needs. It is recommended not only to use adhesives, but also a tray with retentions. For this reason, the inner base of the tray should be equipped with retention grooves and, if need be, the tray border should be pierced with retention holes in order to improve the stability of the functional impression.

When using an adhesive for taking a functional impression, it should only be applied onto the area of the functional edge. For shaping the functional edge, use **Function** (1/2 measuring spoon of each A and B component is sufficient). After approx. 30 seconds' kneading, insert the material into the application syringe enclosed.

Apply **Function** evenly onto the border of the tray, and, if necessary, remodel by hand and insert the tray into the patient's mouth. As the setting period in the mouth is about 60 seconds, all necessary functional movements have to be completed during this period. After that, the material should harden without any influence of dynamic forces. After approx. 2 - 2½ minutes the tray can be removed.

The dorsal edge (Ah-line) shall not be part of the functional edge and should be impressed seperately to optimise the fit of the dentures to make. As an alternative to **Function, Perfect**, a soft putty addition- curing silicone, may be used to shape functional edges. When using this product, the application syringe, which is enclosed with **Function**, has to be ordered separately.

**The base correction impression:**

When taking a base correction impression of the upper jaw and the lower jaw, quite different aspects have to be taken into account. In both cases, it is recommendable to apply, in a first step, adhesive onto the base area up to the functional edge to improve the adhesive strength of the correction material.

**1. The lower jaw with Mandisil:**

**Mandisil** is a special silicone which has specifically been developed for base corrections of the lower jaw. Its consistency and fluidity have been adapted to the compression and jaw conditions of the lower jaw.

Mix 5 cm of each A and B component of **Mandisil** on a mixing pad, apply it onto the functional tray already equipped with functional edges and insert it in the patient's mouth. After a maximum period of 3 minutes, the functional impression, which is now completed, can be removed.

**2. The upper jaw with S4:**

Due to the considerably larger surface of the upper jaw, the compression conditions are completely different to those of the lower jaw. For this reason, it is recommendable to use, in case of well-developped jaw structures, the base correction material **S4** (a fine-flowing correction material without any thixotropic properties). Should the upper jaw has already atrophied, use a mixture of **Mandisil** and **S4** for the base correction of the upper jaw.

**Mandisil** and **S4** for the base correction of the upper jaw.


Physical data:	Labosil blue	Function	Perfect	Mandisil	S 4
Consistency:	putty blue	putty olive green	soft putty green	medium blue green	light blue
Colour:	0,1	0,1	0,1	0,1	0,1
Linear shrinkage (max. %):	90	60	48	48	50
Final hardness (shore A, approx):	30	30	30	30	30
Mixing time (sec.):	30	30	30	30	30
Processing time: (about min:sec from the beginning of mixing):	2:00	2:30	2:30	2:15	2:30
Setting time in mouth (approx. min.;sec.):	< 3:00	3:30	3:00	3:00	3:00
ISO 4823, type:	type 0	type 0	type 0	type 2	type 3
Bisico addition -curing silicones can be stored at least 2 years (as from the date of delivery).					

<b>Delivery Form:</b>					
<b>Labosil blue (90 shore A)</b>	kneadable/putty	<b>Mandisil</b>	mittelließend/medium		
Package with tins a 660 ml (330 ml A-, 330 ml B-comp. + 2 dosing coins)	<b>REF:</b> 01010	Package with tubes a 200 g/150 ml, (100 g A-,100 g B-comp., 2 tube keys)	<b>REF:</b> 01220		
<b>Function</b>	kneadable/putty	<b>S 4</b>	feinfließend/light		
Package with tins a 300 g/208 ml (150 g A-, 150 g B-comp., 10 g Adhesive , 2 dosing coins ,1 application syringe)	<b>REF:</b> 01350	Package with tubes a 200 g/150 ml, (100 g A-, 100 g B-comp., 2 tube keys)	<b>REF:</b> 01260		
<b>Perfect</b>	kneadable/softputty				
Package with tins a 730 g/ 600 ml (365 g A-,365 g B-Comp., 2 dosing coins)	<b>REF:</b> 01310				

Made in Germany!
For dental use only!
Federal law restricts this device to sale by or on the order of a dentist

Rev.-Nr.: 65610-08.2018

Bisico • Bielefelder Dentalsilicone GmbH & Co. KG • Johanneswerkstraße 3 • D-33611 Bielefeld
Tel.: +49 521 8016800 • Fax: +49 521 8016801 • E-Mail: info@bisico.de

## bisico

## Complete denture prosthetics with addition curing silicones

The mixing ratio of **Mandisil** and **S4** is dependant on the progress of atrophy and may be adapted accordingly. The impression is taken in the same way as that of the lower jaw and can be removed after a maximum period of 3 ½ minutes. The vestibule impression with **Perfect**.

After both functional impressions have been closed in correct relation by using hard silicones (e.g. Labosil blue), a “vestibule impression” with **Perfect** can be taken to obtain various static and cosmetic parameters such as the inner lip closing line, the canine tooth position or the sagittal step. Further informations are available on request.

**Recommendations for the preliminary impression:**

In complete denture prosthetics, it is common practice to take a topographical impression with a customised impression tray which is used as preliminary impression. However, such impressions only reach a mediocre degree of precision and are usually extremely over-extended. Besides, in case of an obvious mistake of impression, there is no means of correction. These problems can be avoided by taking the preliminary impression with the silicone **Labosil blue** which allows for making subsequent corrections after it has been removed from the tray. **Labosil blue** is a very hard and particularly stable, kneadable silicone which achieves, by using the Ludwigs technique” rail, exact results, even without the stabilising effect of an impression tray.

To improve the precision, the preliminary impression, taken by using **Labosil blue**, should be completed with functional edges (**Function**) and the base should be corrected (**Mandisil**). Further information about the complete denture prosthetics system developed by Horst Ludwigs and known as “Ludwigs Technique” are available on request. (Also see Bisico catalogue)

**Other information:**

To obtain a model of highest surface quality, wait about 2-3 hours before pouring the model.

**Disinfection:**
Impressions taken by using the material described above can be disinfected in all usual ways.

**Storage:**
Bisico impression silicones shall be used at a room temperatur of about 25°. Larger quantities should be stored at cellar temperature over a longer period of time, avoid temperatures of more than 30° .

**Liability limitations:**

In case of faults in the material or obvious failures, Bisico's liability is limited to replace the faulty material after having tested it. Other claims, particularly those resulting from the impression methods used or from application faults, are excluded from liability. Side effects and interactions:
The use of latex gloves and various haemostatic substances may adversely affect the reaction of the A-silicones.

**Environmental information:**

Bisico impression silicones do not present any ecological problem when they have hardened and small quantities of this material (e. g. material for one impression) may be disposed of with household refuse. Non-cured silicones must be cured (hardened) before being disposed of!

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