Temp S

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Temp S

Special material for fabricating temporary crowns and bridges

Temp S is an especially developed dental Acryl resin on the basis of Isobutylmethacrylat, for the intra-oral fabrication of temporary crowns and bridges. **Temp S** is deliverable in a universal colour, which correspond roughly to the Vita colour A 3 resp. the bio colour 2!

Additionally we can deliver a darkening and a brightening powder. With this the universal colour can be influenced in direction „dark“ or „light“. The colour correction can be made by partly adding to the universal colour and can be changed, according to the part, directly from „very dark“ to „very light“.

Preparation

Before starting the preparation, a situation impression is made. It can be taken with an addition curing silicone, but an impression of alginate or C-silicone is also sufficient.

When using an A-silicone, the temporary can be produced later once more, if necessary.

In the area of missing teeth, a blocking-out of the tooth space as space holder for an interim part is made before the impression. The tooth form will be created during finishing.

In order to stabilize the temporary, it will be cut off especially distinctive interdental areas in the area of the corresponding teeth.

The fabrication of the temporary is made in 4 phases:

Phase 1:

Duration approx. 3 min. Make a saturated mixture by adding powder to the liquid and mix it to avoid colour streaks. Finally the material can be poured into impression taken before the preparation.

Phase 2:

Duration approx. 2:30 min. Now, **Temp S** is a mouldable, adhesive-free mass and will be applied without pressure together with the impression in the patient mouth. Leave it there for 2:30 min.

Phase 3:

Duration approx. 2:30 min. After removal from the mouth, the temporary will be impressed. Cut the remaining edges with a knife and remove the temporary carefully from the impression. Now, the material is elastic and retrieves its initial shape after every deformation.

In order to compensate polymerisation contractions, the temporary will be inserted several times onto the ground and removed again.

Phase 4:

Duration approx. 2 - 3 min. As soon as the material develops reaction heat, it is strong enough, to let it harden completely outside the mouth. Finally, if necessary, it can be easily finished with stone or milling cutter resp. polished.

Reparation of a temporary:

The defect part will be roughened and if necessary some grooves will be cut. Apply new material in the reparation area. After hardening finishing and polishing follow.

Security advices:

The Monomer (liquid) is flammable!

Please close the monomer immediately after use to avoid smell inconvenience.

Please avoid monomer surplus, to avoid irritations and allergic reactions.

The **Temp S** contains Methacrylate.

Temp S liquid is irritating. Please avoid contact with skin, mucous, and eyes.

When contact with skin, wash with water and soap.

When contact with the eyes, wash them immediately and thoroughly with water and contact immediately the doctor.

Avoid spilling the monomer!

Disposal:

Smaller quantities can be disposed with the house refuse, unhardened rests has to be hardened before disposal. (polymerized).

Packaging:

Temp S

Content: 100 g powder, colour universal, 50 ml liquid, 1 mixing cup

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BISICO • Bielefelder Dentsilicone GmbH & Co KG • Johanneswerkstraße 3 • D-33611 Bielefeld