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Progressive impression technique with condensation curing silicones

General information:

Bisico condensation-curing silicones are well-known for their top quality and their universal areas of application. If the material is treated with gloves, it may be advisable to use nitrile gloves as this is less prone to sticking due to a textured surface. The whole spectrum of dental impression techniques is covered by three kneadable and one free-flowing material as well as by the special material **Regidur blau** required for bite registrations. All condensation-curing silicones can be used either with **Universal paste hardener** or with **liquid hardener A**. The dosage instructions are indicated on the packages or on the mixing pads. In general, red **Universal paste hardener** is preferred as it is easy to handle and allows a visual mixing control by developing a homogeneous colour during mixing.

The dosing of the kneadable material is done with a measuring spoon. The free-flowing material is dosed by means of measuring lengths on the mixing pad. A level spoon or a large measuring length of the silicones on the mixing pad are mixed with a small measuring length of the paste hardener.

When using the **liquid hardener A**, a measuring spoon of kneadable silicone and 8 to 10 drops of the **liquid hardener A** have to be mixed. When using the free-flowing material **Exakt-N**, a large measuring length of the basis material has to be mixed with 4 - 5 drops of the **liquid hardener A**.

The dosage indicated in the instructions refers to current reaction times which have proven in dental practices. With a larger or smaller dosage of the hardener, the reaction time can be shortened or prolonged. An over dosage or underdosage of less than 30 % is allowed and does not have any effect on the silicones. A reduction in quality and reaction stability is possible!

Plast / Exakt N:

In general, the well-known double impression technique is applied. The preliminary impression material **Plast** is dosed and mixed according to instructions, then filled in the impression tray and inserted into the patient's mouth. Provided that the hardener is correctly dosed, the intraoral setting time is about 3 minutes. It is recommended to make a test impression to determine the setting time according to the individual dosage habits. After that, the preliminary impression is removed from the patient's mouth. In a next step, it is prepared for the correction impression by cutting off the cervical margins and the interdental spaces as narrow as possible along the non-prepared teeth. Then the impression is filled, by avoiding any bubbles, with the already mixed correction silicone **Exakt N** and immediately inserted into the mouth. The setting time is about 3 ½ minutes depending on the quantity of hardener used. The plaster models should be made 30 - 120 minutes after taking the impression in order to obtain a model of exact dimensions. Impressions taken by using Bisico condensation-curing silicones are easy to handle and can be used for all dental procedures including electroplating.

Lutesil96/Lutesil V:

Lutesil96 is primarily used in dental laboratories, e. g. for blocking out, for making a matrix and similar purposes. The short processing and setting times and the possibility to process it directly at the technician's workplace enables the technician to work in a very economical way. **Lutesil V** is mainly used in the field of total prosthetics, e.g. for producing bite rims and preliminary functional impressions. It can also be used in the Ludwigs technique.

Regidur blau

Regidur blau is used for bite registrations and achieves a final hardness of about 96 shore. On demand, packages including the **liquid hardener A** are also available.

The material is dosed with the measuring spoon, the hardener is added as described on the left and then kneaded for a period of about 30 seconds. After that, a roll is formed by hand which is placed on the lower dental arch. The mouth is closed up to the terminal occlusion position. After 3 minutes' waiting the forming can be started.

The technician can accurately trim the registration by cutting, milling and grinding and is thus able to place plaster models of correct relation in the articulator.

Storage:

The storage temperature should not exceed 25°C!

Contraindications:

In the very rare case of an allergic reaction to a substance used, Bisico condensation-curing silicones should not be used for taking the impression. In this case, there is no problem to use Bisico addition-curing silicones for taking the impression.

Environmental information:

Bisico impression silicones do not present any ecological problem when they have cured (hardened) and small quantities of this material (e.g. the material for one impression) may be disposed of with household refuse. Non-cured silicones must be cured (hardened) before being disposed off!

Disinfection:

Impressions taken by using Bisico condensation-curing silicones can be disinfected with all kinds of disinfectants (e.g. glutaraldehyde solution) available on the market. The manufacturer's instructions for use have to be followed.

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Technical data:	Regidur blau	Lutesil 96/Lutesil V	Plast	Exakt-N
Consistency:	kneadable	kneadable	kneadable	fine-flowing
Colour:	blue	light-grey	white	blue
Final hardness (approx. shore A):	96	96	60	35
Mixing time (seconds):	30	30	30	30
Processing time (seconds):	90	90	90	100
Intraoral setting time (min.:sec.):	3:00	3:00	3:00	3:30
ADA 19, type 2 / ISO 4823, type 1:	type 0	type 0	type 0	type 3

Bisico condensation-curing impression silicones have a shelf life of at least 2 years (from date of delivery)

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bisico® Bielefelder Dentsilicone GmbH @ Co.KG • Johanneswerkstraße 3 • D-33611 Bielefeld