Micro Esthetic flow-viscous /Micro Esthetic Gingiva

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Micro Esthetic flow-viscous /Micro Esthetic Gingiva

Micro Esthetic flow viscous und **Micro Esthetic Gingiva** are Bis-GMA free, light curing, flowable, radiopaque (210%AD), high viscosity Composites. Complies with the guidelines described in DIN EN ISO 4049.

Micro Esthetic flow viscous und **Micro Esthetic Gingiva** are available in syringes.

Colours Micro Esthetik flow-viscous: **A1 / A2 / A3 / A3.5**

Colours Micro Esthetic Gingiva: **G-L (light) / G-M (medium) / G-D (dark) / G-SD (super dark)**

Composition

Monomer matrix: Diurethane dimethacrylate, 1,4-butanediol dimethacrylate

Total filler: 77% by weight (57% by volume) anorganic filler (0.005 – 40 µm)

Indications Micro Esthetic flow viscous:

- Fissure sealing
- Extensive fissure sealing on molars and premolars
- Fillings in Black's class V (cervical caries, root erosions, v-shaped lesions)
- Minimally invasive fillings in Black's classes I, II and III
- Corrections of enamel defects
- Blocking out undercuts
- Small corrections of shape and color of the enamel

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- Restorations on the gingiva
- Covering of exposed tooth neck
- Fillings in Black's class V (cervical caries, root erosions, v-shaped lesions)

Procedure:

Pre-treatment:

Before starting the treatment, clean the hard tooth substance with a fluoride-free polishing paste. Select the correct color for **Micro Esthetik flow-viscous**.

1. Cavity preparation

Prepare the cavity according to the general rules of the adhesive technique. Use a procedure that is gentle on the hard tooth substance. Bevel the incisal edges on anterior teeth. Do not bevel the edges on posterior teeth. Avoid irregular edges. Next, clean the cavity with water spray, remove all residue and dry. It is necessary to dry the site completely. We recommend the use of a kofferdam.

In case of non carious defects the tooth sureface has to be roughening before using **Micro Esthetic Gingiva**.

2. Pulp protection/liner

The use of a liner can be foregone, if an enamel dentin adhesive is used. In case of very deep cavities that come close to the pulp, cover the according surfaces with a calcium hydroxide compound.

3. Design of the approximal contacts

For cavities in approximal areas, insert and fix a transparent matrix.

4. Adhesive system

Etch (z. B. **Best Etch**) and bond (z. B. **Bond LC**) according to the manufacturer's instructions.

For **Micro Esthetic Gingiva** we recommend a selfetching adhesive, e.g. **E-Bond LC**. Please use the adhesive according to the manufacturer's instructions.

5. Application

Fill the composite in thin layers (max. 2 mm) directly into the cavity. Use the disposable curved application tips included with the material. For hygiene reasons, do not use the application tips more than one time. Avoid introducing air bubbles into the material during application with the tip. Ensure that the prepared hard tooth substance is well-coated. Polymerize each layer with a conventional polymerization lamp for 40 sec. A dispersion layer will form on the surface of the material during polymerization. This dispersion layer forms the chemical bond between the layers and must not be touched or removed.

6. Finishing

The composit is ready to be finished and polished directly after polymerization. Suitable tools are finishing diamonds, flexible wheels, silicone polishers and polishing brushes.

Special notes:

- In case of extensive restorations, the surgical light should be moved away from the working area temporarily to avoid premature curing of the composite. Alternatively, the material can be covered with a light-tight foil.
- For hygiene reasons, do not use the curved application tips included with the product more than once. Dispose of the tips directly after use.
- Use a light polymerization system with an emission range of 350 – 500 nm to polymerize the material. The required physical properties are only reached if the polymerization light functions properly. Therefore, it is necessary to check the light intensity regularly according to the manufacturer's instructions

Light intensity for polymerization : 650 mW/cm²

Wave length for polymerization : 350 – 500 nm

Polymerization time: 40 sec.

Hazard statement / precautionary statement: Contains Tetramethylene dimethylacrylate

Warning:

May cause an allergic skin reaction. Wear protective gloves/protective clothing/eye protection/face protection. If skin irritation or rash occurs: Get medical advice/attention.

Contra-indications

If a patient has known hypersensitivities towards a component of this product, we recommend not to use it or to do so only under strict medical supervision. The dentist should consider known interactions and crossreactions of the product with other materials already in the patient's mouth before using the product.

Side-effects:

With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (allergies) or local discomfort, however, cannot be ruled out completely. Should you learn about unwanted side-effects - even if it is doubtful that the side-effect has been caused by our product - please kindly contact us.

In case of exposed dentin in a cavity, ensure a proper protection of the pulp to avoid possible pulp reactions (e. g. line the cavity with a calcium hydroxide compound).

Cross-reactions with other substances:

Phenolic substances (e. g. eugenol) inhibit polymerization. Therefore, these materials (e. g. zinc oxide eugenol cements) must not be used as liners.

Disinfection / Protection from cross-contamination:

Place the functioning syringe with attached application tip into a suitably shaped barrier sheath; pierce end of sheath with cannula, exposing the cannula for use. Using a barrier sheath facilitates cleaning and disinfection of the syringe between patients. After use of sheathed syringe, remove delivery tip and sheath by grasping on the hub of the delivery tip through the sheath, twist and remove tip along with sheath. Discard used tip and sheath in appropriate waste stream. Replace syringe storage cap. Disinfect - After removal of the delivery tip and sheath, disinfect this product using an intermediate level disinfection process (liquid contact) as recommended by the Centers for Disease Control and endorsed by the American Dental Association, Guidelines for Infection Control in Dental Health-Care Settings – 2003 (Vol. 52; No. RR-17), Centers for Disease Control and Prevention (USA).

Storage:

Store at temperatures between 10 °C and 25 °C (50 °F to 77 °F). Avoid direct sunlight. Screw the cap back onto the syringe tightly after each use. Let the material reach room temperature before use. Withdraw the plunger slightly after use to keep the outlet from becoming plugged. Do not use use after the expiration date (see label on the syringe). Only for use in dentistry. Keep out of children's reach. This product has been developed for the specific use illustrated above. Only process as described in these instructions. The manufacturer will not be held liable for any damages that result from improper use or improper processing.

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PT

Delivery Form:

Micro Esthetic flow-viscous, A1 Syringe with 2 g, 5 applicationtips	REF 75075
Micro Esthetic flow-viscous, A2 Syringe with 2 g, 5 applicationtips	REF 75070
Micro Esthetic flow-viscous, A3 Syringe with 2 g, 5 applicationtips	REF 75071
Micro Esthetic flow-viscous, A3.5 Syringe with 2 g, 5 applicationtips	REF 75072
Additional products: Application tips, 1,1 mm x 13 mm bended, 20 pieces	REF 75150
Micro Esthetic Gingiva, G-L (light) Syringe with 2 g, 5 applicationtips	REF 75076
Micro Esthetic Gingiva, G-M (medium) Syringe with 2 g, 5 applicationtips	REF 75077
Micro Esthetic Gingiva, G-D (dark) Syringe with 2 g, 5 applicationtips	REF 75078
Micro Esthetic Gingiva, G-SD (super dark) Syringe with 2 g, 5 applicationtips	REF 75079

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*Vita is a registered trademark of Vita® Zahnfabrik H. Rauter GmbH & Co. KG, Bad Säckingen
* Vita é uma marca registrada da Vita* Zahnfabrik H. Rauter GmbH & Co. KG, Bad Säckingen
bisico® Bielefelder Dentsilsilicone GmbH @ Co.KG • Johanneswerkstraße 3 • D-33611 Bielefeld