Cem DC

Dual Cure Comosite Cement

Cem DC is a dual cure easy flowing microhybrid composite for permanent cementations of inlays, onlays, crowns, bridges and veneers. **Cem DC** is a **strongly radiopaque** and relatively high filled composite of very high strength.

Cem DC is based on methacrylate resin and inorganic filler particles of **0.05-1 µm**. The total filler load is 61 % and the total filler volume 41 %. Cem DC meets the requirements of **DIN EN ISO 4049**.

Preparing the Syringe

Remove the cap of the syringe and throw it away (**do not use it again!**). It is replaced by a special 1:1 mixing cannula. Turn the cannula 90° until it locks in position. The Material is now ready for application.

The working time (23°C (74°F)) in the self cure mode is 2:00 minutes from start of mixing.

Note:

Store used syringe with fixed used mixing cannula. Discard the first 2-3 mm (about the size of a peppercorn) of the extruded material. This has to be done for each new mix.

Cementation of Inlays, Onlays, Crowns, Bridges and Veneers**1. Preparing**

Prepare the luting side areas of the restoration (inlay, onlay, crown, bridge or veneer) with a suitable bonding agent according to the corresponding instructions.

2. Cementing of the indirect restoration**2.1. Cementing of Inlays, Onlays, Crowns and Bridges**

Prepare the tooth surface with a suitable dual cure bonding agent (**E-Bond DC**). For cementing (preparation of the restoration see 1.1.) apply a uniform coating (0,5 mm) of **Cem DC** on the luting side of the inlays, onlays, crowns and bridges and on the primed tooth areas. Seat the restoration under light pressure, and retain it in place until the cement has completely cured. Remove excess material. Light cure all marginal areas of the restoration from each direction for **20 seconds** with a polymerization unit (wavelength range 400-500 nm) with a light intensity of at least 1000mW/cm². The dual cure cement system will auto cure within **4 minutes**.

2.2. Cementing of Veneers

Optional try-in the veneers with Try-in Gel. Clean try-in Gel off the enamel surfaces, using pumice in a rubber cap. Avoid gingival contact to prevent bleeding. Rinse with water and dry with oil-free air. Isolate the teeth to be veneered with interproximal strips to protect adjacent teeth (not being veneered) from the etchant and bonding agent.

Apply an etchant (e.g. **Best Etch**) onto the enamel areas to be veneered. Leave the etchant in place for 15 seconds, and rinse with plenty of water. After drying place a thin coat of suitable dual curing adhesive system (ega mixture of **Bond LC** and **Activator SC**) as described in the corresponding user instruction before application of **Cem DC**.

Remove matrix strips prior to placing veneers on teeth. This step assures complete and passive seating of the veneer, even in multiples.

After placement of the veneer, light cure for **20 seconds** with a polymerization unit (wavelength range 400-500 nm) with a light intensity of at least 1000mW/cm² through tooth structures and the veneer material. If the light through the veneer material is not sufficient or a light cure cannot guaranteed to be sufficient, **Cem DC** selfcures in **4 minutes**. After curing finish and polish margins in the usual manner.

Additional Notes

- Do not use any resin to adjust viscosity of composite restorative material.
- Do neither store the composite material in proximity of eugenol containing products, nor let the composite allow to come into contact with materials containing eugenol. Eugenol can impair the hardening of the composite
- Contact of resin pastes with skin should be avoided, especially by anyone having known resin allergies.

Storage

Do not store above 20 °C (68 °F). Store unopened material in the refrigerator.

Opened cartridges have been used up within 3 months.

Does not use after expiry date.

Warranty

Bisico GmbH warrants this product will be free from defects in material and manufacture. Bisico makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusively remedy and Bisico's sole obligation shall be repair or replacement of the Bisico product.

Limitation of Liability

Except where prohibited by law, Bisico GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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