## Multi Repair Set

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
Rev.: 65820-04.2021

## Multi Repair Set

The combination of the Multi Repair Primer and Multi Repair Bond represents a fast and universally applicable system for intraoral and extraoral repair of fractured ceramic blending without acid cauterization.

### Ceramic Repair System

#### Preparation of Surface Fractures

Smooth the damaged surface and roughen it with a fine grit diamond (< 50 µm) or intraoral air polisher before cleaning off any adhering material with oil free compressed air. Avoid contamination with moisture.

#### Please note

The use of an intraoral air polisher (Al2O3/< 120 µm) further enhances the bond strength!

#### The further repair procedue will result in 3 steps:

##### Step 1: Multi Repair Primer

In combination with Multi Repair Bond, Multi Repair Primer serves as an adhesive agent for ceramic/metal and composites.

Using a brush (white shank), apply a thin coat of Multi Repair primer to the ceramic or any exposed metal immediately after cleansing with oil free air pressure. Allow to air dry for approximately 1 minute. The primer must be applied onto the whole repair surface; we recommend a small overfilling of the surface. A second coat may be applied if necessary. The treated surface must be dry; otherwise adhesion will be reduced considerably.

##### Step 2: Multi Repair Opaker

Multi Repair Opaker serves as a cover for the metal frame. Therefore, its use is necessary only when the fracture extends past the metal frame.

Using a hard, shorthaired brush (red shank), apply a **thin** coat of the Multi Repair Opaker. A poorly covering Multi Repair Opaker coat adversely influences the color effect. At least two coats are necessary to achieve total coverage. The fluidity of the Multi Repair Opaker is optimized through stirring on a mixing block. Special attention must be given to the thin coating in the retention area. Each coat is polymerized separately (see table for polymerization time ranges).

Thick opaque coatings prevent optimal polymerization, therefore weakening the bond. A longer polymerization time does not lead to greater hardening depth.

#### Processing range of the Opaker

1 - 3 minutes, depending on light conditions

##### Step 3: Multi Repair Bond

Multi Repair Bond is not intended for completion of the blending. It merely improves the adhesion bond between the ceramic and the composite, and may only be applied thinly. The actual repair is then carried out with a composite material.

#### Application

Multi Repair Bond is applied onto the fracture surface with a needle and, using a brush, is distributed into a thin coat to cover the entire surface (disposable single use needles included).

Multi Repair Bond is of gel from consistency and can only be distributed with the enclosed disposable brushes (red shank) with red shank correctly and thinly (this brush is also suitable for the application of the Multi Repair Opaker).

Then the paste is hardened with a polymerization light (either halogen or LED tools are suitable). The processing time of Multi Repair Bond ranges between 30 - 60 seconds (depending on light conditions).

#### Polymerization times for Multi Repair Opaker and Multi Repair Bond with different light hardening tools (valid as of: 10/2004) (data not guaranteed)

Polymerization times	Multi Repair Opaker	Multi Repair Bond
<b>Light hardening tool</b>		
Spektra 2000 (Schütz Dental)	180 sec.	180 sec.
Spektra LED (Schütz Dental)	50 sec.	20 sec.
Uni XS (Heraeus Kulzer)	180 sec.	180 sec.
Optilux 501 (Kerr)	60 sec.	40 sec.
CapoLux (Schütz Dental Translux CL)	60 sec.	40 sec.
(Heraeus Kulzer)	60 sec.	40 sec.

#### Light hardening of composite fillings

Apply composite according to manufacturer's instructions, harden with light, process and polish.

To avoid air inclusions, the composite should be modeled strongly and specially on the transitions to the ceramic blending very carefully.

#### Important!

Care must be taken to ensure that neither saliva nor other contaminants come into contact with the surface being treated during the entire repair procedure, since this would considerably reduce the adhesion bond. We recommend the use of a rubber dam in case of large restorations!

Avoid any contact of Multi Repair Primer and Multi Repair Bond with the mucous membranes.

To avoid early hardenings on the products, you should strongly reduce or switch off the operation light of the treatment unit.

#### Storage

Store the material at 5-25°C. Close product immediately after use and protect from direct light.

#### Shelf life

The maximum shelf life is printed on the label of each product. Please do not use after the expiration date.

#### Troubleshooting list

Errors	Cause	Remedy
Blending does not hold after repair	Surface not sufficiently roughened, cleaned and dried Multi Repair Opaker or Multi Repair Bond not entirely hardened faulty frame dimensions the repaired composite was not placed out of occlusion	Rough out frame surfaces and remove all grease layers, moisture and milling dust see following points new dimensioning of frame in areas exposed to severe masticatory loading, trim composite to prevent occlusal contact
Multi Repair Opaker not completely hardened	Multi Repair Opaker applied too thickly Single coats not polymerized separately Insufficient light output by polymerization lamp	apply thinner coat of Multi Repair Opaker separate polymerization of individual coats check light output and eventual change of light source
Multi Repair Bond not completely hardened (paste still green after light hardening)	Multi Repair Bond applied too thickly Insufficient light output by polymerization lamp	apply thinner coat of Multi Repair Bond check light output and eventual change of light source
Multi Repair Bond in the needle appears to be too hard and solid	Material not correctly sealed, paste prepolymerized Material not warmed to room temperature after removal from the refrigerator	Replace cap on needle after each material withdrawal Allow Multi Repair Bond to warm to room temperature before use
Multi Repair Bond in the needle appears too fluid	Material was stored above 25°C for a longer time	Observe the storage temperature, storage in refrigerator if necessary

#### Warranty

Our technical user tips, regardless of whether they are provided verbally, in writing or through practical guidance, stem from our own experiences and can therefore be regarded only as guidelines. Our products undergo continuous development. Therefore, we reserve the right to possible future changes.

#### Indications

Please observe the existing relevant safety data sheets while processing.

#### Ingredients

Primer: Methylmetacrylat

Bond: Activator containing (Campherchinon) liquid of Methacrylacidester

Opaker: Activator containing (Campherchinon) liquid of Methacrylacidester

#### Side effects

With appropriate processing and use of this medical product, undesired side effects are rarely experienced. Immune reactions (such as allergies) or local irritations cannot be ruled out entirely in principle. Please inform us of any undesired side effects you may have experienced - even when in doubt.

#### Adverse signs/interactions

Do not use or use only under strict supervision of the attending physician/dentist if patient suffers from hypersensitivity. In this instance, the composition of our medicinal product is available upon request.

Known cross reactions or interactions of this medicinal product with other materials already in the mouth must be taken into account by the dentist prior to use.

0482

Rev.: 65820-04.2021